



adaptiQ® Sterile Injection Vials

SCHOTT Pharma ISO standard sterile and depyrogenated, ready-to-use vials, supplied in a choice of two packaging formats



adaptiQ® Vials in Cup Nest

SCHOTT Injection Vials supplied sterile nestled in cups within a plastic preform, inside an industry standard tub. Cup Nest packaging offers simple de-nesting and re-nesting and is also ideal for automated processing.



adaptiQ® Vials in a Tray

Sterile vials in a tray are ideal for large production runs using traditional infeed equipment, or for situations where full removal of the vials is required for manual filling.

Cup Nest

NBB (No Blowback)

Size	Adelphi Part No.	Quantity per Tub	MOQ
2R	VC2RNBBBQC	100	15 Tubs
6R	VC6RNBBBQC	48	15 Tubs
10R	VC10RNBBBQC	48	15 Tubs
20R	VC20RNBBBQC	24	15 Tubs
50R	VC50RNBBBQC	16	15 Tubs

Tray

NBB (No Blowback)

Size	Adelphi Part No.	Quantity per Tray	MOQ
2R	VC2RNBBBQT	160	10 Trays
6R	VC6RNBBBQT	88	10 Trays
10R	VC10RNBBBQT	80	8 Trays
20R	VC20RNBBBQT	48	8 Trays
50R	VC50RNBBBQT	28	6 Trays

Unsure which option is best suited to your process?
Contact our friendly sales team and we will be happy to assist you:

sales@adelphi-hp.com